



Disposable HMEFs

一次性使用呼吸过滤器(人工鼻)

新跃医疗旗下品牌

HUAKUN 华坤

Disposable HMEFs

一次性使用呼吸过滤器(人工鼻)

- Multiple filters to satisfy various requirements
- Matched with relative consumables such as anesthesia masks and circuits
- Clinically used for heating, humidifying and filtering during mechanical ventilation of severe and anesthetized patients
- High bacterial viral efficiency >99.99%, insulation and moisturizing
- Single use, 24 hours replacement is recommended
- 提供多种人工鼻可满足临床不同的需求
- 可以与麻醉面罩、管路等相关耗材连接使用
- 临床用于重症及麻醉病人机械通气时气体的加温、湿化和过滤
- 产品具有高效的病毒和细菌过滤率(病毒、细菌过滤率大于99.99%)
- 单人单次使用，建议24小时更换



Model 产品编码	A1(3000)	B2(3006)	F1(3019)
Luer port Luer 接口	√	√	√
Filtration efficiency 过滤效率	>99.99%	>99.99%	>99.99%
Resistance 流速阻力(30L/min)(kPa)	0.09	0.2	0.09
Resistance 流速阻力(60L/min)(kPa)	0.2	0.4	0.19
Compressible volume 可压缩容积(mL)	19	17.5	17.4
Weight 重量(g)	30.6	20.6	32.8
Minimum tidal volume 最小潮气量(mL)	100	50	100

Model 产品编码	F3(3136)	F5(3055)	I1(3321)	G3(3308)
Luer port Luer 接口	√	√	√	√
Filtration efficiency 过滤效率	>99.99%	>99.99%	>99.99%	>99.99%
Resistance 流速阻力(30L/min)(kPa)	0.08	0.19	0.05	0.05
Resistance 流速阻力(60L/min)(kPa)	0.2	0.5	0.11	0.11
Compressible volume 可压缩容积(mL)	19.9	4.1	29	48.6
Weight 重量(g)	30.3	14.1	24.5	25.1
Minimum tidal volume 最小潮气量(mL)	100	50	75	100

宁波华坤医疗器械有限公司

地址：浙江省慈溪市高新技术产业开发区滨江路208号
 电话：0574-6390 5558
 传真：0574-6307 3411
 邮编：315301
 邮箱：huakun@xinhaigroup.com
 官网：www.chinahuakun.cn

Ningbo Huakun Medical Equipment Co., Ltd.

Add: No.208 Binjiang Road, Cixi High-tech Industrial Development
 Zone, Cixi City, Zhejiang, 315301, China
 Tel.: 86-574-6302 9812
 Fax: 86-574-6307 3411
 E-mail: huakun@xinhaigroup.com
 Web: www.chinahuakun.cn